

CoughAssist T70

Prescription

Patient name _____

Date of birth _____

Date _____

Settings

Mode: manual auto

Interface: mask

Cough-Trak: on off

mouthpiece

trach adaptor

Inspiratory pressure (range) _____ cm H₂O

Inspiratory time _____ secs.

Expiratory pressure (range) _____ cm H₂O

Expiratory time _____ secs.

Frequency of treatment _____

Pulse oximetry SpO₂ ≤ _____ % NA

Titrate to achieve inspiratory pressures of 35-40 cm H₂O and expiratory pressures of -35-45 cm H₂O or to achieve an effective cough.

Additional orders

Duration of use

lifetime other _____

Physician information

Name (please print) _____

Signature _____

Telephone _____