

Caregiver competency checklist

Patient name: _____ Caregiver name: _____

Goals	Date goals met	Trainer
The caregiver will be able to:		
1. Understand basic indicator		
2. Review safety Information.		
3. Explain NIOV system components. <ul style="list-style-type: none"> <input type="checkbox"/> Purge tube <input type="checkbox"/> Oxygen regulator <input type="checkbox"/> NIOV <input type="checkbox"/> Pillows interface <input type="checkbox"/> Oxygen Cylinder <input type="checkbox"/> Oxygen supply hose <input type="checkbox"/> Belt Clip <input type="checkbox"/> Battery charger 		
4. Explain the buttons, indicators, and connectors on the front, rear, and side of the NIOV. <ul style="list-style-type: none"> <input type="checkbox"/> Touch Screen <input type="checkbox"/> Power Button <input type="checkbox"/> Speaker <input type="checkbox"/> Charger Connection <input type="checkbox"/> Alarm Silence Button <input type="checkbox"/> Oxygen supply connection <input type="checkbox"/> Activity Buttons <input type="checkbox"/> Power Light <input type="checkbox"/> Backup alarm buzzer <input type="checkbox"/> Breath indicator light <input type="checkbox"/> Interface port 		
5. Attach O2 source and discuss O2 safety.		
6. Assemble and connect: <ul style="list-style-type: none"> <input type="checkbox"/> Purge tube <input type="checkbox"/> Oxygen supply hose <input type="checkbox"/> Pillows interface <input type="checkbox"/> Belt Clip <input type="checkbox"/> Battery charger 		
7. Turn Niov on and off.		



Goals	Date goals met	Trainer
The caregiver will be able to:		
<p>8. Navigate the home and menu screens.</p> <p>Home screen:</p> <p><input type="checkbox"/> Wrench button <input type="checkbox"/> Flip button</p> <p><input type="checkbox"/> Activity icon <input type="checkbox"/> Battery charge icon</p> <p><input type="checkbox"/> Vibration icon <input type="checkbox"/> Time/date icon</p> <p>Menu screen:</p> <p><input type="checkbox"/> Home button <input type="checkbox"/> Settings button</p> <p><input type="checkbox"/> Information button</p>		
<p>9. Locate and explain alarm messages and alarm indicators (audio and visual).</p> <p><input type="checkbox"/> Explain the audio pause feature</p> <p><input type="checkbox"/> Explain what to do when an alarm occurs</p>		
<p>10. Review care, cleaning, and maintenance of NIOV and accessories.</p>		
<p>11. Review emergency procedures.</p>		

Caregiver signature _____ Technician signature _____

